

University of Florida - IFAS
Soil and Water Science Department Core Laboratories
ICP-Mass Spectrometry Laboratory

Mission Statement- The Biogeochemistry of Trace Metals Laboratory (BTML) provides expertise and service for elemental analysis of environmental samples.

Location- 3183 McCarty Hall A

Equipment- NexION 300X inductively coupled plasma mass spectrometer (ICP-MS), Evolution 300 UV-VIS Spectrophotometer.

Procedures and Costs-

1. Elemental analysis:

	UF Accounts	External Pricing
(a) Experimental setup -	\$165	\$220
(b) Analysis of 1 element, per sample -	\$7	\$9
(c) Additional elements, per sample -	\$4	\$5

2. Other procedures:

(a) Digestion/Extraction, per sample -	\$11	\$17
(b) Filtration, per sample -	\$11	\$17

Note: Prices listed may be appropriately discounted for labor or materials provided from external sources. Use of the laboratory facilities by graduate students and post-doctoral associates will be accommodated to the extent possible.

For additional information contact:

Jonathan Judy

Soil and Water Science Department

University of Florida, P.O. Box 110290, Gainesville, FL, 32611

Tel: 352-294-3143; E-mail: jonathan.judy@ufl.edu

Basis for Establishing ICP-Mass Spectrometry Laboratory Analytical Prices for each procedure are established to cover the following costs on a per sample basis for the following:

- (i) Use of expendable materials (tubing, standards, argon, chemicals, tips, gloves etc.).
- (ii) Maintenance and operation costs of equipment required to conduct the procedure.
- (iii) Labor costs for the level of expertise required for the procedure.
- (iv) Overhead and incidental costs, on a time-proportional basis, that have to be paid by the core lab (e.g., electrical service and plumbing; safe chemical waste disposal, etc.).

The following are example cost analysis (UF accounts) for operation of an ICP-MS program:

Analysis of 10 samples, 4 elements: \$165 setup, \$19 per sample: \$355

Analysis of 10 samples, 1 element: \$165 setup, \$7 per sample: \$235

Analysis of 10 samples, 2 elements: \$165 setup, \$11 per sample: \$275

Sample Preparation:

Solutions should be diluted so that the final acid concentration is approximately 1-3%, usually aqua regia or nitric acid (**NOT hydrofluoric**).

Solutions should be filtered so there is no particulate matter.

Minimum volume required is typically 5-10 mL.

Always include a blank digest prepared with the same protocols as your sample.

Let me know if you have any additional questions.